

A new high resolution, eight element Silicon Drift Detector for µPIXE analysis at LIBAF

M. Elfman, L. Ros, N. De La Rosa, P. Kristiansson, E.J.C. Nilsson and J. Pallon. Div. of Nuclear Physics, Lund University, Lund, Sweden



Introduction

- The Lund Ion Beam Analysis Facility (LIBAF) has been upgraded with a high resolution, high count rate 8 x 80 mm² annular Silicon Drift Detector (SDD X-PIPS) [1]) for PIXE analysis. The detector arrangement is a retrofit to the old liquid nitrogen cooled cryostat [2].
- Solid angle increased by tilting the normal of each of the elements towards the



beam line axis. Solid angle covered by the elements increased by about 60% to 265 msr, compared to flat geometry.

- Resolution at 5.9 keV Mn K_{α} line at optimum conditions is 126 129 eV. \bullet
- Modified Mesytec MSCF16F unit [3,4] used as 16 channel spectroscopy amplifier, together with a multihit TDC module allows for the possibility to combine energy signals with timing information, building a software pile-up rejector.



Fig. 1: Principal electronics layout (simplified) of spectroscopy and timing signals. MSCF-16: 16 channel spectroscopy amplifier with remote control, timing filter CFD/LE and a multiplicity discriminator trigger, but no built-in pile up rejection filters. This is instead implemented with a multihit TDC (CAEN V1290N) [3,4], where software cuts are made offline in the TDC-1D or in the 2D plot and data is resorted.

Fig. 3 Sample arrangement



Fig. 2: Efficiency of new SDD (pink) compared with previous HPGe (turquoise). Green is ratio times 0.1





Fig. 4: Left: 2D plot of energy spectrum vs time. Right: PIXE spectrum from thick Ti target. Filter effect can clearly be seen. Green is with and blue without pileup rejection.

Summary

Fig. 5: Spectra from one single detector element and the sum of all eight. All spectra identical and can be summed together directly. A lot of the small peaks are lifted up above the background.

- 8 channel SDD has been successfully implemented
- High count rate keeping the resolution, p/b 1/10000
- High energy resolution (126 eV), directly sum all eight channels
- Short shaping time (0.25 μ s)

References:

[1] CA Canberra Industries, Inc. 800 Research Parkway, Meriden, CT 06450, U.S.A. [2] A. Shariff, et al. Nucl. Instr. Meth. B 219–220 (2004) 494-498 [3] M. Elfman, et al. Nucl. Instr. Meth. B 371 (2016) 148– 85640 Putzbrunn, Germany, (<http://www.mesytec.com>).

Contact: Mikael Elfman/Linus Ros, Div. of Nuclear Physics, Lund University, Lund, Sweden; Mikael.Elfman@nuclear.lu.se, Linus.Ros@nuclear.lu.se [4] Mesytec GmbH & Co. KG, Wernher-von-Braun-Str. 1, 15th International Conference of Nuclear Microprobe Technology and Applications ICNMTA 2016, July 30 - August 5, 2016, Lanzhou, China